




Accesorios para tanque de electroforesis de gel MultiSUB Mini

Cleaver Scientific

Accesorios para tanque de electroforesis de gel MultiSUB Mini



| Code        | Description                  | Packaging |
|-------------|------------------------------|-----------|
| LLG06227144 | Comb 8 sample, 1.5mm thick   | 1 pz.     |
| LLG06227253 | Comb 16 sample, 1mm thick    | 1 pz.     |
| LLG06227254 | Comb 16 sample, 0.75mm thick | 1 pz      |
| LLG07659130 | Comb, 10 sample, 1mm thick   | 1 pz      |
| LLG07671891 | Comb 4 Prep, 2 mm thick      | 1 pz      |
| LLG09584655 | Comb, 8 sample, 1mm thick    | 1 pz.     |
| LLG09584656 | Comb 16 sample, 1.5mm thick  | 1 pz.     |



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

