

**Micropipetas TwiXS Pack Acura<sup>®</sup> manual XS 826****Socorex****Micropipetas TwiXS Pack Acura<sup>®</sup> manual XS 826**

Seis paquetes de pipetas diferentes. Para aplicaciones específicas en el laboratorio de investigación. Relación atractiva calidad/precio.

Un paquete listo para su uso incluye:

- Dos pipetas Acura<sup>®</sup> manual XS
- Un soporte de dos pipetas para repisa
- Muestras de puntas Qualitix<sup>®</sup>
- Certificado de Control de Calidad
- Instrucciones de utilización



Code	Description	Packaging
LLG09411650	Acura® manual 826 XS TwiXS Pack Epack/2	1 pz.
LLG09411651	Acura® manual 826 XS TwiXS Pack Gpack/2	1 pz.
LLG09411652	Acura® manual 826 XS TwiXS Pack Kpack/2	1 pz.
LLG09411653	Acura® manual 826 XS TwiXS Pack Npack/2	1 pz.
LLG09411654	Acura® manual 826 XS TwiXS Pack Hpack/2	1 pz.
LLG09411655	Acura® manual 826 XS TwiXS Pack Qpack/2	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5ª Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

