



**Agujas de disección**

Wironit. CrNi 18/12. Longitud: 50 mm.



Code	Description	Packaging
LLG09160030	Dissecting needles,CrNi-steel,length 50	1 pz.
LLG09160035	Dissecting needles,CrNi-steel,length 50	1 pz.
LLG09160036	Semi lancet, 50 mmpack of 10/10	10 pz.
LLG09160037	Spatula, 50 mmpack of 10/10	10 pz.
LLG09160040	Dissecting needles,CrNi-steel,length 50	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

