



Asas de cultivo, de acero Wironit

Hammacher

Asas de cultivo, de acero Wironit
De acero Wironit al CrNi 18/12.



Code	Description	Packaging
LLG09160055	Loops,CrNi-steel,length 50 mm,diam. 4,0	1 pz.
LLG09160056	Inoculationsloops 2,0mm, Wironit 1/4452	1 pz.
LLG09160057	Loops,CrNi-steel,length 50 mm,diam. 1,5	1 pz.
LLG09160058	Loops,CrNi-steel,length 50 mm,diam. 2,5	1 pz.
LLG09160061	Loops,CrNi-steel,length 60 mm,diam. 1,5	1 pz.
LLG09160062	Loops,CrNi-steel,length 60 mm,diam. 2,5	1 pz.
LLG09160063	Loops,CrNi-steel,length 60 mm,diam. 4,0	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

