




## Accesorios para el tanque de electroforesis de gel MultiSUB Midi



| Code        | Description                   | Packaging |
|-------------|-------------------------------|-----------|
| LLG07657241 | Comb 25 samples, 1,5 mm thick | 1 pz      |
| LLG09584657 | Comb 8 sample, 1mm thick      | 1 pz.     |
| LLG09584658 | Comb 10 sample MC, 1mm thick  | 1 pz.     |
| LLG09584659 | Comb 12 sample, 1mm thick     | 1 pz.     |
| LLG09584660 | Comb 16 sample, 1mm thick     | 1 pz.     |
| LLG09584662 | Comb 25 sample, 1mm thick     | 1 pz.     |
| LLG09584811 | Comb 20 sample MC, 1 mm thick | 1 pz.     |



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5º Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

