

**Tubos de ensayo LLG, vidrio Fiolax[®]****LLG Labware****Tubos de ensayo LLG, vidrio Fiolax[®]**

Con reborde y fondo redondo. Resistencia a las altas temperaturas y a los productos químicos.



Code	Description	Packaging
LLG06234026	LLG-Test tubes, FIOLAX, 10x100mmwit/100	100 pz
LLG07190001	LLG-Test tubes, FIOLAX, 8x70 mmwith/100	100 pz
LLG07190011	LLG-Test tubes, FIOLAX, 12x100 mmwi/100	100 pz
LLG07190021	LLG-Test tubes, FIOLAX, 16x160 mmwi/100	100 pz
LLG07190023	LLG-Test tubes, FIOLAX, 18x180 mmwi/100	100 pz
LLG07190028	LLG-Test tubes, FIOLAX, 20x180 mmwi/100	100 pz
LLG07606374	LLG-Test tubes, FIOLAX, 14x130 mmwi/100	100 pz
LLG07606375	LLG-Test tubes, FIOLAX, 16x130 mmwi/100	100 pz
LLG07613795	LLG-Test tubes, FIOLAX, 12x75 mmwit/100	100 pz



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

