

**Tubos de ensayo LLG, vidrio Fiolax<sup>®</sup>****LLG Labware****Tubos de ensayo LLG, vidrio Fiolax<sup>®</sup>**

Con reborde y fondo redondo. Resistencia a las altas temperaturas y a los productos químicos.



Code	Description	Packaging
LLG06234026	LLG-Test tubes, FIOLAX, 10x100mmwit/100	100 pz
LLG07190001	LLG-Test tubes, FIOLAX, 8x70 mmwith/100	100 pz
LLG07190011	LLG-Test tubes, FIOLAX, 12x100 mmwi/100	100 pz
LLG07190021	LLG-Test tubes, FIOLAX, 16x160 mmwi/100	100 pz
LLG07190023	LLG-Test tubes, FIOLAX, 18x180 mmwi/100	100 pz
LLG07190028	LLG-Test tubes, FIOLAX, 20x180 mmwi/100	100 pz
LLG07606374	LLG-Test tubes, FIOLAX, 14x130 mmwi/100	100 pz
LLG07606375	LLG-Test tubes, FIOLAX, 16x130 mmwi/100	100 pz
LLG07613795	LLG-Test tubes, FIOLAX, 12x75 mmwit/100	100 pz



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

