



Accesorios para tanque de electroforesis en gel MultiSUB Choice

Accesorios para tanque de electroforesis en gel MultiSUB Choice



| Code | Description | Packaging |
|-------------|----------------------------------|-----------|
| LLG06231475 | UV TRAY 15X15 CM | 1 PZ |
| LLG06233618 | Comb 20 samples, 1,5 mm thick | 1 pz. |
| LLG06234451 | Comb 20 samples, 1 mm thick | 1 pz. |
| LLG06238277 | Comb 35 samples, 1 mm thick | 1 pz. |
| LLG06258162 | Comb 16 samples MC, 1.5 mm thick | 1 pz |
| LLG09584665 | Comb 18 samples MC, 1 mm thick | 1 pz. |
| LLG09584666 | Comb 30 samples MC, 1,5 mm thick | 1 pz. |
| LLG09584667 | Comb 18 samples MC, 1,5 mm thick | 1 pz. |



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

