



Accesorios para tanque de electroforesis de gel MultiSUB Maxi

Accesorios para tanque de electroforesis de gel MultiSUB Maxi



Code	Description	Packaging
LLG06251650	MultiSUB Maxi 20 x 20cm UV Tray	1 pz
LLG06269204	Comb, 25 sample, 1,5mm thick	1 pezzo
LLG09584672	Comb 20 samples MC, 1 mm thick	1 pz.
LLG09584673	Comb, 30 sample, 1mm thick	1 pz.
LLG09584675	Comb, 50 sample, 1mm thick	1 pz.
LLG09584676	Comb, 20 sample MC, 1,5mm thick	1 pz.
LLG09584677	Comb, 30 sample, 1,5mm thick	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

