



Accesorios para tanque de electroforesis de gel MultiSUB Maxi

Accesorios para tanque de electroforesis de gel MultiSUB Maxi



Code	Description	Packaging
LLG06251650	MultiSUB Maxi 20 x 20cm UV Tray	1 pz
LLG06269204	Comb, 25 sample, 1,5mm thick	1 pezzo
LLG09584672	Comb 20 samples MC, 1 mm thick	1 pz.
LLG09584673	Comb, 30 sample, 1mm thick	1 pz.
LLG09584675	Comb, 50 sample, 1mm thick	1 pz.
LLG09584676	Comb, 20 sample MC, 1,5mm thick	1 pz.
LLG09584677	Comb, 30 sample, 1,5mm thick	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

