



Accesorios para el tanque de electroforesis de gel VS20 Wave Maxi

Accesorios para el tanque de electroforesis de gel VS20 Wave Maxi



Code	Description	Packaging
LLG06232288	Notched glass plates 200 x 200 mm2 mm/2	1 pz.
LLG06232289	Glass plate 200 x 200 mmplain, with b/2	1 pz.
LLG06232630	Comb, 18 samples MC1 mm thick	1 pz.
LLG06252009	Glass plate 200x200mm, notchedwith 1m/2	2 pz
LLG09584737	Plain glass plates 200x200mmwith 1,5m/2	1 pz.
LLG09584738	Comb with 24 samples1mm thick	1 pz.
LLG09584740	Comb with 48 samples1mm thick	1 pz.
LLG09584741	Comb with 18 samples MC1.5mm thick	1 pz.
LLG09584742	Comb with 24 samples1.5mm thick	1 pz.



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

