



Accesorios para el tanque de electroforesis de gel VS20 Wave Maxi

Accesorios para el tanque de electroforesis de gel VS20 Wave Maxi



Code	Description	Packaging
LLG06232288	Notched glass plates 200 x 200 mm 2 mm/2	1 pz.
LLG06232289	Glass plate 200 x 200 mm plain, with b/2	1 pz.
LLG06232630	Comb, 18 samples MC 1 mm thick	1 pz.
LLG06252009	Glass plate 200x200mm, notched with 1m/2	2 pz
LLG09584737	Plain glass plates 200x200mm with 1,5m/2	1 pz.
LLG09584738	Comb with 24 samples 1mm thick	1 pz.
LLG09584740	Comb with 48 samples 1mm thick	1 pz.
LLG09584741	Comb with 18 samples MC 1.5mm thick	1 pz.
LLG09584742	Comb with 24 samples 1.5mm thick	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

