



Termómetro de propósito general LLG, relleno rojo

LLG Labware

Termómetro de propósito general LLG, relleno rojo
Graduación resistente a ácidos y bases.



Code	Description	Packaging
LLG06206594	Thermometer, red filling, -20/+110°C	1 pz
LLG09236720	LLG-Thermometer, general purpose -35...	1 pz
LLG09236721	LLG-Thermometer, general purpose -10...	1 pz
LLG09236722	LLG-Thermometer, general purpose -10...	1 pz
LLG09236723	LLG-Thermometer, general purpose -20...	1 pz
LLG09236724	LLG-Thermometer, general purpose -10...	1 pz
LLG09236725	LLG-Thermometer, general purpose -10...	1 pz
LLG09236726	LLG-Thermometer, general purpose -10...	1 pz
LLG09236727	LLG-Thermometer, general purpose -10...	1 pz
LLG09236730	LLG-Thermometer, general purpose -10...	1 pz.
LLG09236731	LLG Thermometer, general purpose -20...	1 pz



CARLO ERBA Reagents SA
Calle Filadors 35,
6ª Planta 5ª Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

