



**Termómetro de propósito general LLG, relleno rojo**

**LLG Labware**

**Termómetro de propósito general LLG, relleno rojo**  
Graduación resistente a ácidos y bases.



Code	Description	Packaging
LLG06206594	Thermometer, red filling, -20/+110°C	1 pz
LLG09236720	LLG-Thermometer, general purpose -35...	1 pz
LLG09236721	LLG-Thermometer, general purpose -10...	1 pz
LLG09236722	LLG-Thermometer, general purpose -10...	1 pz
LLG09236723	LLG-Thermometer, general purpose -20...	1 pz
LLG09236724	LLG-Thermometer, general purpose -10...	1 pz
LLG09236725	LLG-Thermometer, general purpose -10...	1 pz
LLG09236726	LLG-Thermometer, general purpose -10...	1 pz
LLG09236727	LLG-Thermometer, general purpose -10...	1 pz
LLG09236730	LLG-Thermometer, general purpose -10...	1 pz.
LLG09236731	LLG Thermometer, general purpose -20...	1 pz



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6ª Planta 5ª Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

