



**Termómetro de precisión LLG, calibrable, relleno rojo**

**LLG Labware**

**Termómetro de precisión LLG, calibrable, relleno rojo**

Capilar: recubierto de amarillo, redondo.



Code	Description	Packaging
LLG09236783	LLG-Precision-Laboratory Thermometer -1	1 pz
LLG09236784	LLG-Precision-Laboratory Thermometer -1	1 pz
LLG09236785	LLG-Precision-Laboratory Thermometer -1	1 pz
LLG09236786	LLG-Precision-Laboratory Thermometer -1	1 pz
LLG09236787	LLG-Precision-Laboratory Thermometer 0.	1 pz
LLG09236788	LLG-Precision-Laboratory Thermometer -1	1 pz
LLG09236800	LLG-Precision-Cold-Laboratory Thermomet	1 pz
LLG09236801	LLG-Precision-Cold-Laboratory Thermomet	1 pz



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

