



Termómetro de precisión LLG, calibrable, relleno rojo

LLG Labware

Termómetro de precisión LLG, calibrable, relleno rojo

Capilar: recubierto de amarillo, redondo.



Code	Description	Packaging
LLG09236783	LLG-Precision-Laboratory Thermometer -1	1 pz
LLG09236784	LLG-Precision-Laboratory Thermometer -1	1 pz
LLG09236785	LLG-Precision-Laboratory Thermometer -1	1 pz
LLG09236786	LLG-Precision-Laboratory Thermometer -1	1 pz
LLG09236787	LLG-Precision-Laboratory Thermometer 0.	1 pz
LLG09236788	LLG-Precision-Laboratory Thermometer -1	1 pz
LLG09236800	LLG-Precision-Cold-Laboratory Thermomet	1 pz
LLG09236801	LLG-Precision-Cold-Laboratory Thermomet	1 pz



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

