



**Termómetro de precisión LLG, calibrable, relleno rojo**

**LLG Labware**

**Termómetro de precisión LLG, calibrable, relleno rojo**

Capilar: recubierto de amarillo, redondo.



| Code        | Description                             | Packaging |
|-------------|---|-----------|
| LLG09236783 | LLG-Precision-Laboratory Thermometer -1 | 1 pz      |
| LLG09236784 | LLG-Precision-Laboratory Thermometer -1 | 1 pz      |
| LLG09236785 | LLG-Precision-Laboratory Thermometer -1 | 1 pz      |
| LLG09236786 | LLG-Precision-Laboratory Thermometer -1 | 1 pz      |
| LLG09236787 | LLG-Precision-Laboratory Thermometer 0. | 1 pz      |
| LLG09236788 | LLG-Precision-Laboratory Thermometer -1 | 1 pz      |
| LLG09236800 | LLG-Precision-Cold-Laboratory Thermomet | 1 pz      |
| LLG09236801 | LLG-Precision-Cold-Laboratory Thermomet | 1 pz      |



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

