



## Tubos de ensayo LLG, vidrio sodocálcico

### Tubos de ensayo LLG, vidrio sodocálcico

Tubos de ensayo con fondo redondo. Resistencia al agua: clase HGB 3 conforme a DIN ISO 719, tipo III conforme a USP.



Code	Description	Packaging
LLG06201034	Test tubes, Ø 8mm, H 70mm, w/beaded rim	100 pz
LLG06267806	LLG-Test tubes 100 x 30 mmAR glass,/100	100 pz
LLG06800382	LLG-Test tubes 130 x 16 mmAR glass,/100	100 pz
LLG06800620	LLG-Test tubes 130 x 14 mmAR glass,/100	100 pz
LLG06800621	LLG-Test tubes 150 x 25 mmAR glass, /50	50 pz
LLG06801073	LLG-Test tubes 100 x 10 mmAR glass,/100	100 pz
LLG07600324	LLG-Test tubes 200 x 30 mmAR glass, /50	50 pz
LLG07600405	LLG-Test tubes 160 x 16 mmAR glass,/100	100 pz
LLG07607198	LLG-Test tubes 100 x 12 mmAR glass,/100	100 pz
LLG07607778	LLG-Test tubes 180 x 20 mmAR glass,/100	100 pz
LLG07651871	LLG-Test tubes 120 x 16 mmAR glass,/100	100 pz.
LLG07651872	LLG-Test tubes 100 x 15 mmAR glass,/100	100 pz.
LLG07652552	LLG-Test tubes 130 x 18 mmAR glass,/100	100 pz.
LLG09400030	Test tubes,soda-glass,130 x 14 mm,pack	100 pz.
LLG09400040	Test tubes,soda-glass,160 x 16 mm,/100	100 pz.
LLG09400050	Tube à essais fond rond, en verre s/100	100 pz.





**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5º Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

