



Tapones de rosca para frascos de cuello estrecho, PP 28

Tapones de rosca para frascos de cuello estrecho, PP 28

Tapones de PP (a prueba de fugas) 28, con boca interior de aprox. 19,4 mm.



| Code | Description | Packaging |
|-------------|--|-----------|
| LLG06272536 | Screw cap DIN 28, blackPP 500 pz | 500 pz |
| LLG06272537 | Closure DIN 28, PPwhite, with cone seali | 500 pz |
| LLG09070211 | Closure DIN 28, PPwhite, with joint/500 | 500 pz. |
| LLG09070213 | Originality closure DIN 28, PPwhite/500 | 500 pz. |



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

