



Tapones de rosca para frascos de cuello estrecho, PP 28

Tapones de rosca para frascos de cuello estrecho, PP 28

Tapones de PP (a prueba de fugas) 28, con boca interior de aprox. 19,4 mm.



Code	Description	Packaging
LLG06272536	Screw cap DIN 28, blackPP 500 pz	500 pz
LLG06272537	Closure DIN 28, PPwhite, with cone seali	500 pz
LLG09070211	Closure DIN 28, PPwhite, with joint/500	500 pz.
LLG09070213	Originality closure DIN 28, PPwhite/500	500 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

