



Caja plegable clax, 46 litros

Caja plegable clax, 46 litros

Robusta caja plegable de plástico de alta calidad para el uso con el carro de laboratorio clax Mobil.



Code	Description	Packaging
LLG09148099	Clax folding box 46l	1 pz.



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

