

**Tapones de tubo de ensayo LLG, PE****LLG Labware****Tapones de tubo de ensayo LLG, PE**

Hechos de polietileno de alta calidad. Ideales para el sellado a prueba de fugas de tubos de ensayo.



Code	Description	Packaging
LLG09231270	LLG-Test tube stoppers, neutralfor /100	100 pz
LLG09231271	LLG-Test tube stoppers, bluefor 12//100	100 pz
LLG09231272	LLG-Test tube stoppers, yellowfor 1/100	100 pz
LLG09231273	LLG-Test tube stoppers, redfor 12/1/100	100 pz
LLG09231274	LLG-Test tube stoppers, greenfor 12/13m	100 pz
LLG09231275	LLG-Test tubes stoppers, neutralfor/100	100 pz
LLG09231276	LLG-Test tube stoppers, bluefor 15//100	100 pz
LLG09231277	LLG-Test tube stoppers, yellowfor 1/100	100 pz
LLG09231278	LLG-test tube stoppers, redfor 15/1/100	100 pz
LLG09231279	LLG-Test tube stoppers, greenfor 15/100	100 pz



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

