



## Tapones de tubo de ensayo LLG, PE

LLG Labware

### Tapones de tubo de ensayo LLG, PE

Hechos de polietileno de alta calidad. Ideales para el sellado a prueba de fugas de tubos de ensayo.



Code	Description	Packaging
LLG09231270	LLG-Test tube stoppers, neutralfor /100	100 pz
LLG09231271	LLG-Test tube stoppers, bluefor 12//100	100 pz
LLG09231272	LLG-Test tube stoppers, yellowfor 1/100	100 pz
LLG09231273	LLG-Test tube stoppers, redfor 12/1/100	100 pz
LLG09231274	LLG-Test tube stoppers, greenfor 12/13m	100 pz
LLG09231275	LLG-Test tubes stoppers, neutralfor/100	100 pz
LLG09231276	LLG-Test tube stoppers, bluefor 15//100	100 pz
LLG09231277	LLG-Test tube stoppers, yellowfor 1/100	100 pz
LLG09231278	LLG-test tube stoppers, redfor 15/1/100	100 pz
LLG09231279	LLG-Test tube stoppers, greenfor 15/100	100 pz



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

