


**Homogenizadores, cuerpo redondo, vidrio****GPE****Homogenizadores, cuerpo redondo, vidrio**

Homogenizadores completos de vidrio esterilizable, con un contenido de 5 hasta 55 ml. Versiones de pistilo con rosca para accionamiento por motor o bien huecos en su interior para su uso con refrigerantes, por ejemplo, hielo seco. Versiones de mortero como vidrio transparente o vidrio esmerilado para homogeneizar más eficazmente tejidos resistentes. Espacio entre pistilo y mortero 0,152 hasta 0,254 mm



| Code | Description | Packaging |
|-------------|---|-----------|
| LLG09164660 | Homogenizer 5mlround, clear glass,pestl | 1 pz |
| LLG09164661 | Homogenizer 10mlround, clear glass,pest | 1 pz |
| LLG09164662 | Homogenizer 15mlround, clear glass,pest | 1 pz |
| LLG09164663 | Homogenizer 30mlround, clear glass,pest | 1 pz |
| LLG09164664 | Homogenizer 55mlround, clear glass,pest | 1 pz |
| LLG09164665 | Homogenizer 5mlround, clear glass | 1 pz |
| LLG09164666 | Homogeniyer 10mlround, clear glass | 1 pz |
| LLG09164667 | Homogenizer 15mlround, clear glass | 1 pz |
| LLG09164668 | Homogenizer 30mlround, clear glass | 1 pz |
| LLG09164669 | Homogenizer 55mlround, clear glass | 1 pz |
| LLG09164670 | Homogenizer 5mlround, grounded,pestle h | 1 pz |
| LLG09164671 | Homogenizer 10mlround, grounded,pestle | 1 pz |
| LLG09164672 | Homogenizer 15mlround, grounded,pestle | 1 pz |
| LLG09164673 | Homogenizer 30mlround, grounded,pestle | 1 pz |
| LLG09164674 | Homogenizer 55mlround, grounded,pestle | 1 pz |
| LLG09164675 | Homogenizer 5mlround, grounded | 1 pz |
| LLG09164676 | Homogenizer 10mlround, grounded | 1 pz |
| LLG09164677 | Homogenizer 15mlround, grounded | 1 pz |
| LLG09164678 | Homogenizer 30mlround, grounded | 1 pz |
| LLG09164679 | Homogenizer 55mlround, grounded | 1 pz |





CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

