


**Homogenizadores, cuerpo redondo, vidrio****GPE****Homogenizadores, cuerpo redondo, vidrio**

Homogenizadores completos de vidrio esterilizable, con un contenido de 5 hasta 55 ml. Versiones de pistilo con rosca para accionamiento por motor o bien huecos en su interior para su uso con refrigerantes, por ejemplo, hielo seco. Versiones de mortero como vidrio transparente o vidrio esmerilado para homogeneizar más eficazmente tejidos resistentes. Espacio entre pistilo y mortero 0,152 hasta 0,254 mm



Code	Description	Packaging
LLG09164660	Homogenizer 5mlround, clear glass,pestl	1 pz
LLG09164661	Homogenizer 10mlround, clear glass,pest	1 pz
LLG09164662	Homogenizer 15mlround, clear glass,pest	1 pz
LLG09164663	Homogenizer 30mlround, clear glass,pest	1 pz
LLG09164664	Homogenizer 55mlround, clear glass,pest	1 pz
LLG09164665	Homogenizer 5mlround, clear glass	1 pz
LLG09164666	Homogeniyer 10mlround, clear glass	1 pz
LLG09164667	Homogenizer 15mlround, clear glass	1 pz
LLG09164668	Homogenizer 30mlround, clear glass	1 pz
LLG09164669	Homogenizer 55mlround, clear glass	1 pz
LLG09164670	Homogenizer 5mlround, grounded,pestle h	1 pz
LLG09164671	Homogenizer 10mlround, grounded,pestle	1 pz
LLG09164672	Homogenizer 15mlround, grounded,pestle	1 pz
LLG09164673	Homogenizer 30mlround, grounded,pestle	1 pz
LLG09164674	Homogenizer 55mlround, grounded,pestle	1 pz
LLG09164675	Homogenizer 5mlround, grounded	1 pz
LLG09164676	Homogenizer 10mlround, grounded	1 pz
LLG09164677	Homogenizer 15mlround, grounded	1 pz
LLG09164678	Homogenizer 30mlround, grounded	1 pz
LLG09164679	Homogenizer 55mlround, grounded	1 pz





CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

