

**Microscopio educativo, BA81****Motic****Microscopio educativo, BA81**

- Tubo monocular (BA81A-MS corded)/tubo binocular (BA81B-MS w/100X corded), inclinación 45°, capaz de girar 360°
- Ocular WF 10X/18 mm
- Portaobjetivos revolver cuadruple, dirigido hacia atras
- Portaobjetivos acromático EA 4x, 10x, 40x (muelle)
- Portaobjetivos acromático EA 4x,10x,40x (muelle) y 100x (muelle/aceite)(BA81B-MS w/100x corded)
- accionamiento fino, grueso coaxial, ajustable
- Carro en cruz coaxial, integrado
- N-A . 1,25 Abbe Kondensar enfocable
- Diafragma iris con soporte de filtro
- Iluminación LED regulable 20 mA, 3,5V, 70 mV
- Alimentación 220-240V, enchufe VDE (CE)
- Carcasa protectora del polvo



Code	Description	Packaging
LLG09727033	Educational microscope BA81B-MS w/100x	1 pz.



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5º Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

