




## Estantes para congeladores verticales

### Estantes para congeladores verticales

Estantes de acero inoxidable, adecuados para todas las marcas de congeladores verticales. Para cajas de hasta 136 mm x 136 mm.



| Code        | Description                              | Packaging |
|-------------|------------------------------------------|-----------|
| LLG09698950 | Rack, 4x3 boxes 50 mm 221 x 419 x 139/1  | 1 pz      |
| LLG09698954 | Rack, 2x3 boxes 75 mm for upright freez  | 1 pz.     |
| LLG09698959 | Rack, 6x4 boxes 50 mm for upright freez  | 1 pz.     |
| LLG09698960 | Rack, 4x4 boxes 78 mm 556 x 140 x 326 mm | 1 pz.     |
| LLG09698961 | Rack, 5x4 boxes 78 mm 556 x 140 x 407 mm | 1 pz.     |
| LLG09698969 | Rack, 4x4 boxes 50 mm 558 x 139 x 221 mm | 1 pz.     |
| LLG09698970 | Rack, 3x4 boxes 50 mm 562 x 139 x 168 mm | 1 pz.     |
| LLG09698971 | Rack, 5x4 boxes 50 mm 558 x 139 x 276 mm | 1 pz.     |
| LLG09698972 | Rack, 7x4 boxes 50 mm 558 x 139 x 386 mm | 1 pz.     |
| LLG09698973 | Rack, 4x5 boxes 50 mm 688 x 139 x 221 mm | 1 pz.     |
| LLG09698974 | Rack, 5x3 boxes 50 mm 419 x 139 x 276 mm | 1 pz.     |
| LLG09698975 | Rack, 3x4 boxes 75 mm 558 x 139 x 240 mm | 1 pz.     |
| LLG09698976 | Rack, 3x3 boxes 75 mm 419 x 139 x 240 mm | 1 pz.     |



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6ª Planta 5ª Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*



