




Estantes para congeladores verticales

Estantes para congeladores verticales

Estantes de acero inoxidable, adecuados para todas las marcas de congeladores verticales. Para cajas de hasta 136 mm x 136 mm.



Code	Description	Packaging
LLG09698950	Rack, 4x3 boxes 50 mm 221 x 419 x 139/1	1 pz
LLG09698954	Rack, 2x3 boxes 75 mm for upright freez	1 pz.
LLG09698959	Rack, 6x4 boxes 50 mm for upright freez	1 pz.
LLG09698960	Rack, 4x4 boxes 78 mm 556 x 140 x 326 mm	1 pz.
LLG09698961	Rack, 5x4 boxes 78 mm 556 x 140 x 407 mm	1 pz.
LLG09698969	Rack, 4x4 boxes 50 mm 558 x 139 x 221 mm	1 pz.
LLG09698970	Rack, 3x4 boxes 50 mm 562 x 139 x 168 mm	1 pz.
LLG09698971	Rack, 5x4 boxes 50 mm 558 x 139 x 276 mm	1 pz.
LLG09698972	Rack, 7x4 boxes 50 mm 558 x 139 x 386 mm	1 pz.
LLG09698973	Rack, 4x5 boxes 50 mm 688 x 139 x 221 mm	1 pz.
LLG09698974	Rack, 5x3 boxes 50 mm 419 x 139 x 276 mm	1 pz.
LLG09698975	Rack, 3x4 boxes 75 mm 558 x 139 x 240 mm	1 pz.
LLG09698976	Rack, 3x3 boxes 75 mm 419 x 139 x 240 mm	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*



