

**Soporte para embudos, acero cromado****Isolab****Soporte para embudos, acero cromado**

Altura ajustable 140 - 230 mm. Diámetro interior del anillo superior: 45 mm.

Diámetro exterior del anillo inferior: 126 mm.



Code	Description	Packaging
LLG04008510	Funnel holder, stainless steel for 140 -	1 pz.



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5º Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

