

**Soporte para embudos, acero cromado****Isolab****Soporte para embudos, acero cromado**

Altura ajustable 140 - 230 mm. Diámetro interior del anillo superior: 45 mm. Diámetro exterior del anillo inferior: 126 mm.



Code	Description	Packaging
LLG04008510	Funnel holder, stainless steel for 140 -	1 pz.



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

