



Pinzas para cubreobjetos LLG, con cierre automático, acero inoxidable

LLG Labware

Pinzas para cubreobjetos LLG, con cierre automático, acero inoxidable
Acero inoxidable 4301.



Code	Description	Packaging
LLG09160380	LLG-Forceps for cover slips, self-locki	1 pz
LLG09160381	LLG-Forceps for cover slips, self-locki	1 pz



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

