



Cubetas para tinción

Cubetas para tinción

Vidrio sodocálcico. Con tapa encajable.



Code	Description	Packaging
LLG09163273	Staining cylinders,soda glass,with cove	1 pz.
LLG09163275	Staining cylinders,soda glass,with cove	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

