



## Cubetas para tinción

### Cubetas para tinción

Tipo Hellendahl.

Para 16 portaobjetos de 76 x 26 mm.



Code	Description	Packaging
LLG09163277	Staining troughs,glass,Hellendahl,for/1	1 pz.
LLG09163278	Staining troughs,glass,Hellendahl,for/1	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

