



Cubetas para tinción

Cubetas para tinción

Tipo Hellendahl.

Para 16 portaobjetos de 76 x 26 mm.



Code	Description	Packaging
LLG09163277	Staining troughs,glass,Hellendahl,for/1	1 pz.
LLG09163278	Staining troughs,glass,Hellendahl,for/1	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

