



**Pinzas LLG, acero 18/10**

**LLG Labware**

**Pinzas LLG, acero 18/10**  
Punta roma-recta.



| Code        | Description                             | Packaging |
|-------------|-----------------------------------------|-----------|
| LLG04008471 | LLG-Forceps 105 mm, blunt/straightgener | 1 pz      |
| LLG04008472 | LLG-Forceps 115 mm, blunt/straightgener | 1 pz      |
| LLG04008473 | LLG-Forceps 130 mm, blunt/straightgener | 1 pz      |
| LLG04008474 | LLG-Forceps 160 mm, blunt/straightgener | 1 pz      |
| LLG04008475 | LLG-Forceps 200 mm, blunt/straightgener | 1 pz      |
| LLG04008476 | LLG-Forceps 250 mm, blunt/straightgener | 1 pz      |
| LLG04008477 | LLG-Forceps 300 mm, blunt/straightgener | 1 pz      |
| LLG04008478 | LLG-Forceps 115 mm, blunt/straightgener | 1 pz      |
| LLG04008479 | LLG-Forceps 130 mm, blunt/straightgener | 1 pz      |
| LLG04008480 | LLG-Forceps 145 mm, blunt/straightgener | 1 pz      |
| LLG04008481 | LLG-Forceps 160 mm, blunt/straightgener | 1 pz      |
| LLG06255826 | LLG-Forceps 145 mm, blunt/straightgener | 1 pz      |



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5ª Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

