



Pinzas o abrazaderas para agitador MaxQ (excepto MaxQ 8000)

Thermo Scientific

Pinzas o abrazaderas para agitador MaxQ (excepto MaxQ 8000)



Code	Description	Packaging
LLG04009332	Clamp for 25 ml Erlenmeyer flask	1 pz.
LLG04009335	Clamp for 2000 ml Erlenmeyer flask	1 pz.
LLG04009336	Clamp for 4000 ml Erlenmeyer flask	1 pz.
LLG04009337	Clamp for 6000 ml Erlenmeyer flask	1 pz.
LLG04009338	Clamp for 2500 ml Culture flasklow form	1 pz.
LLG04009339	Clamp for 2800 ml Fernbach flask	1 pz.
LLG04009340	Clamp for Microplate-/Deep-well plate	1 pz.
LLG04010229	Clamps for Erlenmeyer Flask 10ml	1 pz.
LLG04010230	Clamps for Erlenmeyer Flask 50ml	1 pz.
LLG04010231	Clamps for Erlenmeyer Flask 250ml	1 pz.
LLG04010232	Clamps for Erlenmeyer Flask 500ml	1 pz.
LLG04010233	Clamps for Erlenmeyer Flask 1000ml	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

