




Racks o gradillas para agitador MaxQ

Thermo Scientific

Racks o gradillas para agitador MaxQ



Code	Description	Packaging
LLG04009341	Angle tube rack, red6x6 boring, 10-13 m	1 pz.
LLG04009342	Angle tube rack, orange6x12 boring, 14-	1 pz.
LLG04009343	Angle tube rack, red6x6 boring, 14-16 m	1 pz.
LLG04009344	Angle tube rack, white4x10 boring, 17-2	1 pz.
LLG04009345	Angle tube rack, white4x5 boring, 17-20	1 pz.
LLG04009346	Angle tube rack, blue4x10 boring, 21-25	1 pz.
LLG04009347	Angle tube rack, blue4x5 boring, 21-25	1 pz.
LLG04009348	Angle tube rack, green3x8 boring, 26-30	1 pz.
LLG04009349	Angle tube rack, green3x3 boring, 26-30	1 pz.
LLG04009350	Rack for 1.5 ml micro tubes8x12 boring,	1 pz.
LLG04009351	Rack for 1.5 ml micro tubes4x6 boring,	1 pz.
LLG04010234	Test Tube Racks 10-13 mm, 6x12, red	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

