



Cabezales de bombas para PLP 380, PLP 2000 y PLP 6000

Behr

Cabezales de bombas para PLP 380, PLP 2000 y PLP 6000



Code	Description	Packaging
LLG09700435	Pump head PPH 1032200 ml/min, for PLP 3	1 pz.
LLG09700436	Pump head PPH 3031600 ml/min, for PLP 3	1 pz.
LLG09700437	Pump head PPH 506148 ml/min, for PLP 38	1 pz.
LLG09700438	Pump head PPH 506248 ml/min, for PLP 38	1 pz.
LLG09700442	Pump head PPH 6036000 ml/min, for PLP 6	1 pz.
LLG09700443	Pump head PPH 70421800 ml/min, for PLP	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

