



Guía de luz flexible para fuentes de luz fría SCHOTT

Guía de luz flexible para fuentes de luz fría SCHOTT



Code	Description	Packaging
LLG09705168	Flexible light guide, 1-arm M-PVC, Ø /1	1 pz.
LLG09705225	Flexible light guide, 1-arm,length 600	1 pz.
LLG09705226	Flexible light guide, 2-arm,length 1000	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

