



Abrazaderas para la fuente de luz fría KL 300 LED y el punto de luz EasyLED

Abrazaderas para la fuente de luz fría KL 300 LED y el punto de luz EasyLED



Code	Description	Packaging
LLG06237972	Mounting clamp for Ø 32 mm for KL 200/1	1 pz.
LLG07910430	Mounting clamp for Ø 25 mm column,for/1	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

