



Tamices para molino seco POLYMIX[®] PX-MFC 90 D

Tamices para molino seco POLYMIX[®] PX-MFC 90 D



Code	Description	Packaging
LLG04013109	Sieve 0.5 mm (for PX-MFC 90 D and PX-/1	1 pz.
LLG04013110	Sieve 0.8 mm (for PX-MFC 90 D and PX-/1	1 pz.
LLG04013111	Sieve 1.5 mm (for PX-MFC 90 D and PX-/1	1 pz.
LLG04013112	Sieve 2.0 mm (for PX-MFC 90 D and PX-/1	1 pz.
LLG04013113	Sieve 3.0 mm (for PX-MFC 90 D and PX-/1	1 pz.
LLG04013114	Sieve 4.0 mm (for PX-MFC 90 D and PX-/1	1 pz.
LLG04013115	Sieve 5.0 mm (for PX-MFC 90 D and PX-/1	1 pz.
LLG04013116	Sieve 6.0 mm (for PX-MFC 90 D and PX-/1	1 pz.
LLG06261186	Sieve 1.0 mm (for PX-MFC 90 D and PX-/1	1 pz.
LLG07980251	Sieve 0.2 mm (for PX-MFC 90 D and PX-/1	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

