



Tarro Cilindrico para Preparaciones, DURAN<sup>®</sup>

DURAN Group

**Tarro Cilindrico para Preparaciones, DURAN<sup>®</sup>**

Tapa esmerilada con asa tipo bola.



Code	Description	Packaging
LLG07021920	JAR SPECIMEN WITH KNOBBED LID 110X253MM	1 pz.
LLG09173416	Specimen jars,DURAN®,with knob lid,80x1	1 pz.
LLG09173417	Specimen jars,DURAN®,with knob lid,120x	1 pz.
LLG09173424	Specimen jars, DURAN, with knob lid, 15	1 pz.
LLG09173426	Specimen jars,DURAN®,with knob lid,100x	1 pz.
LLG09173428	Specimen jars,DURAN®,with knob lid,150x	1 pz.
LLG09173438	Specimen jars,DURAN®,with knob lid,100x	1 pz.
LLG09173439	Specimen jars,DURAN®,with knob lid,150x	1 pz.
LLG09173450	Specimen jars,DURAN®,with knob lid,150x	1 pz.
LLG09173457	Specimen jars,DURAN®,with knob lid,120x	1 pz.
LLG09173459	Specimen jars, 300x250mm, DURAN, with k	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

