



Accesorios para fotómetros Lovibond[®]

Accesorios para fotómetros Lovibond[®]



Code	Description	Packaging
LLG06226018	Verification Standard KIT AL450ZERO/430	1 pz.
LLG09699244	Funnel, plastic,with handle, white	1 pz.
LLG09699246	Verification Standard KIT AL400ZERO/430	1 pz.
LLG09699264	Adapter for round vials 16mm	1 pz.
LLG09699266	Stirring rod 13cm, plastic	1 pz.
LLG09699267	Ceaning brush 10cm	1 pz.
LLG09699273	Mains adapter TR15RA150 Level Vfor Mult	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

