



## Accesorios para fotómetros Lovibond<sup>®</sup>

Accesorios para fotómetros Lovibond<sup>®</sup>



Code	Description	Packaging
LLG06226018	Verification Standard KIT AL450ZERO/430	1 pz.
LLG09699244	Funnel, plastic,with handle, white	1 pz.
LLG09699246	Verification Standard KIT AL400ZERO/430	1 pz.
LLG09699264	Adapter for round vials 16mm	1 pz.
LLG09699266	Stirring rod 13cm, plastic	1 pz.
LLG09699267	Ceaning brush 10cm	1 pz.
LLG09699273	Mains adapter TR15RA150 Level Vfor Mult	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

