

**Electronic multichannel microliter pipettes Eppendorf Xplorer® plus, variable**

The electronic pipettes Eppendorf Xplorer plus were specially designed for high professional standards to provide optimal support for you in your work, with a new intuitive operating concept and design based on the proven Eppendorf PhysioCare® Concept.

- Autodispense
- Manual pipetting, multi-aspirate, dilution, sequential pipetting, sequential dispensing, reverse pipetting
- Programming of multiple modes in a sequence
- Setting of fixed volumes
- Password protection
- History function
- Service interval warning



Code	Description	Packaging
LLG09283187	Electronic multichannel microliter pipettes Eppendorf Xplorer® plus, variable, 0.5 - 10 µl	1 pz
LLG09283188	Electronic multichannel microliter pipettes Eppendorf Xplorer® plus, variable, 0.5 - 10 µl	1 pz
LLG09283189	Electronic multichannel microliter pipettes Eppendorf Xplorer® plus, variable, 5 - 100 µl	1 pz
LLG09283190	Electronic multichannel microliter pipettes Eppendorf Xplorer® plus, variable, 5 - 100 µl	1 pz
LLG09283191	Electronic multichannel microliter pipettes Eppendorf Xplorer® plus, variable, 15 - 300 µl	1 pz
LLG09283192	Electronic multichannel microliter pipettes Eppendorf Xplorer® plus, variable, 15 - 300 µl	1 pz
LLG09283193	Electronic multichannel microliter pipettes Eppendorf Xplorer® plus, variable, 50 - 1200 µl	1 pz
LLG09283194	Electronic multichannel microliter pipettes Eppendorf Xplorer® plus, variable, 50 - 1200 µl	1 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

