



Estantes perforados para Incubadoras Heratherm®;

Thermo Elect.LED GmbH
(Kendro)

Estantes perforados para Incubadoras Heratherm™

Acero inoxidable. Incluye 2 soportes de estantes.



Code	Description	Packaging
LLG09534170	Perforated stainless steel shelf for IG	1 pz



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

