

**Puntas de pipeta Eppendorf Combitips advanced[®]****Eppendorf****Puntas de pipeta Eppendorf Combitips advanced[®]**

Para líquidos viscosos de hasta 200 mPa*s y alta presión de vapor, por ejemplo etanol y acetona.

-Codificadas por color con un anillo de color

-Escala de volumen continuo

-Reconocimiento automático de puntas en Multipipette[®] 4780, Multipipette[®] plus, Multipipette[®] stream y Xstream



Code	Description	Packaging
LLG09283135	Combitips® advanced 0,1 mlcolour co/100	100 pz.
LLG09283136	Combitips® advanced 0,2 mlcolour co/100	100 pz.
LLG09283137	Combitips® advanced 0,5 mlcolour co/100	100 pz.
LLG09283138	Combitips® advanced 1,0 mlcolour co/100	100 pz.
LLG09283139	Combitips® advanced 2,5 mlcolour co/100	100 pz.
LLG09283140	Combitips® advanced 5,0 mlcolour co/100	100 pz.
LLG09283141	Combitips® advanced 10,0 mlcolour c/100	100 pz.
LLG09283142	Combitips® advanced 25,0 mlcolour c/100	100 pz.
LLG09283143	Combitips® advanced 50,0 mlcolour c/100	100 pz.
LLG09283144	Combitips® advanced 0,1 mlPCR clean/100	100 pz.
LLG09283145	Combitips® advanced 0,2 mlPCR clean/100	100 pz.
LLG09283146	Combitips® advanced 0,5 mlPCR clean/100	100 pz.
LLG09283147	Combitips® advanced 1,0 mlPCR clean/100	100 pz.
LLG09283148	Combitips® advanced 2,5 mlPCR clean/100	100 pz.
LLG09283149	Combitips® advanced 5,0 mlPCR clean/100	100 pz.
LLG09283150	Combitips® advanced 10,0 mlPCR clea/100	100 pz.
LLG09283151	Combitips® advanced 25,0 mlPCR clea/100	100 pz.
LLG09283152	Combitips® advanced 50,0 mlPCR clea/100	100 pz.
LLG09283153	Combitips® advanced 0,1 mlbiopur, c/100	100 pz.
LLG09283154	Combitips® advanced 0,2 mlbiopur, c/100	100 pz.



Code	Description	Packaging
LLG09283155	Combitips® advanced 0,5 mlbiopur, c/100	100 pz.
LLG09283156	Combitips® advanced 1,0 mlbiopur, c/100	100 pz.
LLG09283157	Combitips® advanced 2,5 mlbiopur, c/100	100 pz.
LLG09283158	Combitips® advanced 5,0 mlbiopur, c/100	100 pz.
LLG09283159	Combitips® advanced 10,0 mlbiopur, /100	100 pz.
LLG09283160	Combitips® advanced 25,0 mlbiopur, /100	100 pz.
LLG09283161	Combitips® advanced 50,0 mlbiopur, /100	100 pz.
LLG09283167	Assortment pack Combitips® advanced(1 p	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

