

**Puntas de pipeta Eppendorf Combitips advanced ®****Eppendorf SE****Puntas de pipeta Eppendorf Combitips advanced ®**

- Para líquidos viscosos de hasta 200 mPa* s y alta presión de vapor, por ejemplo etanol y acetona.
- Codificación de colores con 1 anillo
- Escala de volumen continuo
- Reconocimiento automático de puntas en Multipette® 4780, Multipette® plus, Multipette® stream y Xstream



Code	Description	Packaging
LLG04668754	Combitips® advanced 0,1 ml colour code w	400 pz.
LLG04668755	Combitips® advanced 0,2 ml colour code l	400 pz.
LLG04668756	Combitips® advanced 0,5 ml colour code v	400 pz.
LLG04668757	Combitips® advanced 1,0 ml colour code y	400 pz.
LLG04668758	Combitips® advanced 2,5 ml colour code g	400 pz.
LLG04668759	Combitips® advanced 5,0 ml colour code b	400 pz.
LLG04668760	Combitips® advanced 10,0 ml colour code	400 pz.
LLG04668761	Combitips® advanced 25,0 ml colour code	400 pz.
LLG04668762	Combitips® advanced 50,0 ml colour code	400 pz.
LLG04668763	Combitips® advanced 0,1 ml biopur, colou	400 pz.
LLG04668764	Combitips® advanced 0,2 ml biopur, colou	400 pz.
LLG04668765	Combitips® advanced 0,5 ml biopur, colou	400 pz.
LLG04668766	Combitips® advanced 1,0 ml biopur, colou	400 pz.
LLG04668767	Combitips® advanced 2,5 ml biopur, colou	400 pz.
LLG04668768	Combitips® advanced 5,0 ml biopur, colou	400 pz.
LLG04668769	Combitips® advanced 10,0 ml biopur, colo	400 pz.
LLG04668770	Combitips® advanced 25,0 ml biopur, colo	400 pz.
LLG04668771	Combitips® advanced 50,0 ml biopur, colo	400 pz.





CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

