



Accesorios para medidor de oxígeno SD 315 Oxi

Accesorios para medidor de oxígeno SD 315 Oxi



Code	Description	Packaging
LLG06272364	Electrode, 1,5m cable for SD310 OXI/1	1 pz.
LLG06272365	Electrode, 10m cable for SD310 OXI/1	1 pz.
LLG06272366	Electrode, 30m cable for SD310 OXI/1	1 pz.
LLG06272367	Service kit for oxygen sensor SD310 O/1	1 pz



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

