



Accesorios para LLG-uniPERISTALTICPUMP 1 y 3

LLG Labware

Accesorios para LLG-uniPERISTALTICPUMP 1 y 3



| Code        | Description                              | Packaging |
|-------------|--|-----------|
| LLG06263542 | Pumphead for LLG-uniPERIPUMP 1 with Phar | 1 pz      |
| LLG06263543 | Pumphead for LLG-uniPERIPUMP 1 with Phar | 1 pz      |
| LLG06263544 | Pumphead for LLG-uniPERIPUMP 1 with Phar | 1 pz      |
| LLG06263546 | PharMed® tubing ID 1 mm with connector ( | 1 pz.     |
| LLG06263547 | PharMed® tubing ID 2 mm with connector ( | 1 pz.     |
| LLG06263548 | PharMed® tubing ID 3 mm with connector ( | 1 pz.     |
| LLG06263552 | PharMed® tubing (ID:3,1mm, OD: 6,3mm)    | 1 pz.     |
| LLG06263554 | PharMed® tubing (ID:6,4mm, OD: 9,6mm) wi | 1 pz.     |



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

